

DEP ES

DEP ES System
can promote
absorption of active
cosmetic substances
applied on skin.

Safe, easy to use and
quick, ensures astonishing
results lowering skin
surface with brightness
and consistency.

*The new
technology for
aesthetic skin
treatments.*



DEP ES

DEP ES System device can promote the proper absorption of cosmetic products applied on skin surface both sole or, best, after a microdermabrasion treatment with Aesthi Peel microabrasion system.

DEP ES System has a digital Delivery Display which enables the customer to rule electrical pulses and vibrations intensity and to visualize skin impedance values, thanks to the integrated impedancimeter, allowing to adapt treatments on patients' skin peculiarity.

DEP ES System guarantees, during each treatment, a pleasant and relaxing massage thanks to the vibrant hand-piece, useful both for body and face cure.

DEP ES System jointly with **Aesthi Peel** microdermabrader produces a synergic effect of both peeling and transdermal delivery technologies with astonishing results.

What makes DEP ES System unique?

- Augmentation of active cosmetic compounds absorption all over the body thanks to its ergonomic tips.
- Relaxing massage.
- Constant skin impedance monitoring to obtain optimized treatments for each patient.
- No alterations in delivered drugs pH.

DEP ES System
Is born from
Mattioli Engineering
experience in aesthetical
medicine field with
Dermoelectroporation®
but uses low pulses
intensities limiting their
action on skin surface to
make it more bright and
young.



MATTIOLI

ENGINEERING
AESTHETIC
DIVISION

MATTIOLI ENGINEERING ITALIA S.P.A

Via Petrarca 80-82
50041 Calenzano (Firenze) - Italy
Tel. 055 882247
Fax 055 8874879
info@mattioliengineering.com

N.Verde 800.389.385

www.mattioliengineering.com